---

**DENTSPLY**  
DETREY

## DeTrey® Conditioner 36

**Phosphoric acid gel for the etching of tooth enamel and conditioning of dentine**

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Dansk

# DeTrey® Conditioner 36

## Phosphoric acid gel for the etching of tooth enamel and conditioning of dentine

**DeTrey® Conditioner 36** is a blue-tinted gel containing 36% phosphoric acid for the etching of tooth enamel and conditioning of dentine prior to bonding with resin-based materials.

### COMPOSITION

- Phosphoric acid
- Highly dispersed silicon dioxide
- Detergent
- Pigment
- Water

### INDICATIONS

- Enamel etch technique in connection with composite and compomer restorations
- Total etch technique in connection with composite and compomer restorations.

### CONTRAINDICATIONS

Application to dentine close to the pulp (less than 1 mm).

### WARNINGS

1. DeTrey Conditioner 36 contains 36% phosphoric acid. Causes burns. Avoid contact with oral tissues, eyes and skin. If accidental contact occurs, flush affected area with generous amounts of water. In case of contact with the eyes, immediately rinse with plenty of water and seek medical attention.
2. DeTrey Conditioner 36 gel should extrude easily: **DO NOT USE EXCESSIVE FORCE**. Replace original cap of DeTrey Conditioner 36 tightly after each use to avoid evaporation. Discard needle after use, as needles may clog if gel is allowed to dry inside.
3. Make sure the needle tip is securely fastened to the syringe prior to application of DeTrey Conditioner 36.

### INTERACTIONS WITH DENTAL MATERIALS

Some liners and bases may be etched by phosphoric acid. In general, this does not impair their barrier function.

### STEP-BY-STEP INSTRUCTIONS

#### 1. Cleaning

Clean unprepared enamel with a rubber cup and pumice or a cleaning paste such as Nupro® Prophylaxis Pastes. Wash thoroughly with water spray and air-dry.

Clean freshly cut enamel and dentine with water spray and then air-dry.

## 2. Acid Conditioning of Enamel and Dentine (Total Etch Technique)<sup>1</sup>

It is recommended to follow the Total Etch Technique described below:

Remove cap from the end of the syringe. Attach disposable needle tip to end of the syringe. **Turn tip clockwise to assure that it is fully seated.**

Gently extrude gel directly to surface/s to be treated. Needle tip may be bent for easy access. Start to apply DeTrey Conditioner 36 gel (36% phosphoric acid) to the cavity surfaces at the enamel margins. For best results, condition enamel for at least 15 seconds and dentine for 15 seconds or less.

Discard needle after use as needle may clog if gel is allowed to dry inside. **Replace original cap.**

**CAUTION: Replace cap tightly after each use to avoid evaporation. Gel should extrude easily: DO NOT USE EXCESSIVE FORCE.**

## 3. Rinsing and Drying

Remove gel with aspirator tube and /or vigorous water spray and rinse conditioned areas thoroughly for at least 15 seconds.

Further treatment of the conditioned area depends on the adhesive material which is applied following the conditioning procedure. Consult respective directions for use.

Once the surfaces have been properly treated, they must be kept uncontaminated. If salivary contamination occurs, thoroughly clean with vigorous water spray, dry, and repeat conditioning procedure of enamel for 5 seconds only. Rinse as described above.

Proceed with adhesive procedure.

## STORAGE

Replace cap immediately after use.

Keep out of sunlight. Not to be stored at temperatures exceeding 24 °C.

## BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

### Manufacturer:

**DENTSPLY DeTrey GmbH**

**De-Trey-Str. 1**

**78467 Konstanz**

**GERMANY**

**Phone +49 (0) 75 31 5 83-0**

### Distributor:

**DENTSPLY Limited**

**Hamm Moor Lane**

**Addlestone, Weybridge**

**Surrey KT15 2SE**

**Phone (0 19 32) 85 34 22**

© DENTSPLY DeTrey 2000-08-22

<sup>1</sup> Alternatively to the Total Etch Technique, the conventional Enamel Etch Technique can be followed. In this case, the enamel margins only are treated with DeTrey Conditioner 36 for at least 15 seconds. Then rinse and dry as described below.

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