---

**DENTSPLY**  
DETRAY

Dyract®

**Compomer Restorative System**

**Kompomer-Restorationssystem**

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# Dyract<sup>®</sup>

## Compomer Restorative System

**Dyract<sup>®</sup>** Restorative System combines the major benefits of glass-ionomer cements – adhesion to tooth substance without the need of acid etching, fluoride release and biocompatibility – with the easy handling of a light-curing composite.

**Dyract** is a single-component radiopaque restorative material in 10 tooth shades (Vita<sup>®1</sup> A2, A3, A3.5, A4, B1, B3, O-B3, C2, C3, C4). The material is pre-dosed in Compules<sup>®</sup> Tips.

**Caution:** For dental use only.

### COMPOSITION

#### Dyract

- UDMA resin
- TCB resin
- Strontium-fluoro-silicate glass
- Strontium fluoride
- Photo initiators
- Stabilisers

#### Prime&Bond<sup>®</sup> NT

- Di- and trimethacrylate resins
- Functionalised amorphous silica
- PENTA (dipentaerythritol penta acrylate monophosphate)
- Photoinitiators
- Stabilisers
- Cetylamine hydrofluoride
- Acetone

### INDICATIONS

1. Restoration of class V cavities and cervical lesions.
2. Restoration of class III cavities.
3. Restoration of all classes of cavities in deciduous teeth.
4. Temporary restoration of class I and II cavities of permanent posterior teeth.
5. Restoration of small class I and II cavities (width less than 1/3 of intercuspal distance) of permanent posterior teeth where the restoration has no direct occlusal contact.
6. Base under class I and II composite restorations (open and/or closed sandwich).
7. Core build-up for metal or PFM inlays, onlays, crowns.
8. Luting of ceramic/composite inlays, onlays, crowns (ultrasonic technique).

Dyract expands slightly over time. The potential benefits of this phenomenon (e.g. compensation for polymerisation shrinkage) and possible risks (e. g. stress build-up) have not yet been fully investigated. Therefore, until the results of on-going studies confirm the suitability of the product for other indications, Dyract should only be used for the above listed restorative procedures.

---

<sup>1</sup> Vita is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co.

## CONTRAINDICATIONS

- Direct or indirect pulp capping.
- Core build-up for full ceramic crowns.
- Use in patients with a known allergy to dimethacrylate resins.

## WARNINGS

1. Prime&Bond NT contains methacrylates which may be irritating to skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of soap and water. The product may cause sensitisation by skin contact in susceptible persons. If skin sensitisation occurs discontinue use.
2. Avoid contact of Prime&Bond NT with mucous membranes. After accidental contact, wash and rinse with plenty of water.
3. Prime&Bond NT contains acetone. Acetone is highly flammable. Keep away from sources of ignition – no smoking. Do not breathe vapour. Take precautionary measures against static discharges.
4. Dyract light-cured compomer restorative contains methacrylates which may be irritating to skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of soap and water. The product may cause sensitisation by skin contact in susceptible persons. If skin sensitisation occurs discontinue use.

## PRECAUTIONS

Avoid Prime&Bond NT saturating gingival retraction cord. If Prime&Bond NT soaks into the cord, it may set hard and bond the cord to the underlying tooth surface making removal difficult.

## INTERACTIONS WITH DENTAL MATERIALS

Eugenol containing dental materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.

If H<sub>2</sub>O<sub>2</sub> has been used to clean the cavity, proper rinsing is essential. Higher concentration H<sub>2</sub>O<sub>2</sub> may interfere with the setting of polymerisable material and should not be used prior to the application of Prime&Bond NT.

Prolonged and intensive contact with acetone-containing products may lead to minute dissolution of the outermost surface of calcium hydroxide materials. This has no detrimental effect on the adhesion to the cavity walls.

## ADVERSE REACTIONS

The following adverse reaction has been associated with the use of acetone solutions and acrylate monomers:

- Reversible inflammatory changes of the oral mucosa after accidental contact.

## STEP-BY-STEP INSTRUCTIONS

### 1. Shade Selection

Shade selection should be made prior to the restorative procedure whilst the teeth are hydrated. Remove any extraneous plaque or surface stain. Use the Dyract shade guide provided which contains samples of original Dyract restorative. The colour coding dot on the shade guide matches the coloured cap on the Compules Tip.

Alternatively, a Vita Lumin® Vacuum shade guide may be used. The Dyract shade corresponds to the central part of the respective Vita tooth.

## 2. Cavity Preparation

In all classes of cavity this may be kept to the minimum required for caries removal. Preparation of mechanical retention is not needed except when the restoration is liable to be subjected to high displacement forces.

## 3. Cleaning

Cavity cleanliness is paramount for the development of adhesion.

In cases where no cavity preparation has been made, clean the tooth surface with a rubber cup and pumice or a prophylaxis-paste like Nupro®. Preparing a fresh surface with a finishing bur will significantly increase bond strength to enamel.

Wash surface thoroughly with air/water spray.

Remove rinsing water by blowing gently with an air syringe or blot-dry with a cotton pellet.

Do not desiccate the dentine structure.

## 4. Pulp Protection

For direct or indirect pulp-capping protect the dentine close to the pulp (< 1 mm) with a hard-setting calcium hydroxide liner (e. g. Dycal®), leaving the remaining cavity surface free for bonding with Prime&Bond NT.

## 5. Application of Prime&Bond NT

One layer of Prime&Bond NT is applied:

1. Dispense Prime&Bond NT directly onto a fresh Applicator Tip<sup>2</sup> or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish<sup>2</sup> or standard dappen dish.
2. Immediately apply ample amounts of Prime&Bond NT to thoroughly wet all tooth surfaces. This surface should be saturated which may necessitate additional application of Prime&Bond NT.
3. Leave the surface undisturbed for 20 seconds.
4. Remove solvent by blowing with air from a dental syringe for at least 5 seconds. Surface should have a uniform, glossy appearance. If not, repeat steps 2 to 4.
5. Light-cure for a minimum of 10 seconds<sup>3</sup>. Ensure uniform exposure of all cavity surfaces.
6. Immediately place Dyract compomer over the cured Prime&Bond NT.

## 6. Placement of Dyract

Insert Compules Tip into the notched opening of the applicator gun barrel.

Dispense Dyract directly into the cavity preparation. In deep cavities, incremental placement and curing (in 3 mm layers or less) is recommended to minimise polymerisation shrinkage.

## 7. Curing

Cure each increment separately with a VLC dental polymerisation unit for at least 40 seconds. The tip of the light guide should be held as close as possible to the restoration during curing.

Important: Be sure to expose each area of the entire restoration to the curing light.

Additionally, the restoration should be cured through lingual or buccal enamel walls.

---

<sup>2</sup> DENTSPLY Applicator Dish and Applicator Tips are available from your dental dealer.

<sup>3</sup> When using a high-performance unit such as the curing lights manufactured by DENTSPLY, a curing time of 10 seconds is sufficient. For curing lights with an output lower than 300 mW/cm<sup>2</sup>, a curing time of 20 seconds should be observed.

## **8. Finishing**

Begin finishing immediately after curing. Gross excess material may be removed with fluted finishing burs or diamonds. Finishing is best achieved by using Enhance™ Finishing and Polishing Discs and interproximal finishing and polishing strips. A high final lustre can be obtained by applying Prisma® Gloss™ and Prisma Gloss Extrafine Polishing Pastes.

### **MAINTENANCE OF APPLICATOR GUN**

The applicator gun is sterilisable by autoclave or cold sterilisation solution following the manufacturers' instructions.

It is recommended that the applicator gun be disassembled for assured sterilisation. Partially close the applicator gun and place thumb under the rear portion of the hinge. Push upward and lift hinge separating the applicator gun, exposing the plunger. Remove residual compomer with a soft paper tissue and a suitable solvent (70% alcohol). To reassemble, insert plunger into applicator gun barrel, press components together and snap hinge mechanism in place.

### **STORAGE**

The Prime&Bond NT bottle and Dyract Compules Tips should be tightly closed immediately after use.

Keep out of sunlight.

Not to be stored at temperatures exceeding 25 °C.

Keep Prime&Bond NT in a well ventilated place.

Humidity can adversely affect the properties of unsealed Compules Tips. Therefore keep Compules Tips sealed in their blister pack until use. Under normal ambient conditions, unsealed Compules Tips stay usable for about 4 weeks.

### **BATCH NUMBER AND EXPIRY DATE**

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

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© DENTSPLY DeTrey 1999-04-22
**DENTSPLY DeTrey GmbH**