DENTSPLY

Aquasil Soft Putty – Regular Set
Aquasil Hard Putty – Fast Set

**Very High Viscosity Precision Impression Materials
based on quadrafunctional hydrophilic siloxanes**

Directions for Use	2
Gebrauchsanweisung	6
Mode d'emploi	10
Istruzioni per l'uso	14
Instrucciones de uso	18
Bruksanvisning	22
Instrukcja użytkownika	26
Brugsanvisning	30
Οδηγίες χρήσης	34

Aquasil Soft Putty – Regular Set

Aquasil Hard Putty – Fast Set

Very High Viscosity Precision Impression Materials based on quadrafunctional hydrophilic siloxanes

Aquasil Impression Materials are suitable for all precision impression techniques where excellent hydrophilic properties, dimensional accuracy, high tear strength, and resistance to permanent deformation are needed.

Aquasil Impression Materials comply with the requirements of ISO 4823 for dental elastomeric impression materials.

INGREDIENTS

Base:

Silica Crystalline
 Vinylidimethicone
 Zeolith
 Highly dispersed silica (AEROSIL®¹)
 Paraffine oil
 Polydimethyl-hydrogensiloxane
 Green pigment
 (Cosmetic Green, C.I. 77288)
 Yellow pigment (Arc Yellow A16 N)

Catalyst:

Silica Crystalline
 Vinylidimethicone
 Zeolith
 Highly dispersed silica (AEROSIL¹)
 Paraffine oil
 Platinum complex + Pt-metal
 Calcium carbonate
 1,3-Divinyltetramethyldisiloxane
 Titanium dioxide

INDICATION FOR USE

Aquasil Soft Putty is a regular setting very high viscosity type impression material with optimum flexibility for easy removal from undercuts and is the ideal tray material for the **double mix technique**. It may also be used as primary impression material for the putty/wash technique.

Aquasil Hard Putty is a fast setting, very high viscosity type impression material with ideal mechanical properties as tray material for the **putty/wash technique**. It may also be used as tray material for the double mix technique for quadrant impressions.

CONTRAINDICATIONS

Do not use Aquasil Hard Putty in combination with full arch impression trays for the double mix technique. The low flexibility of the set material may interfere with tray removal.

PRECAUTIONS AND WARNINGS

Remove unset Aquasil Impression Material from clothing with a suitable cleaning solvent.

¹ AEROSIL is a registered trademark of Degussa AG.

INTERACTIONS WITH OTHER MATERIALS

1. Astringents based on aluminium salts may interfere with the setting reaction of vinyl polysiloxanes. Please consult instructions for use of astringent to be used.
2. Aquasil Impression Materials should not be mixed with, or used in conjunction with polyether materials, condensation cured silicones or polysulfide.
3. 'Sulphur-containing' polymeric gloves may interfere with the setting reaction of Aquasil Impression Materials.
 - Do not wear such gloves when kneading Aquasil Soft Putty.
 - Do not touch retraction cords. Handle cords with college pliers.
 - Do not touch prepared tooth.

PREPARATION AND APPLICATION

Tray Preparation

Silfix Adhesive should be used in all techniques with all types of impression trays. Brush a thin layer of Silfix onto tray and allow to dry for 3 minutes.

Preparation of Impression Area

Clean Impression area with air/water spray and remove excess water.

Mixing and filling the tray

Take equal amounts of catalyst and base using the colour coded measuring spoons and knead with fingertips until colour of mix is uniform. Immediately fill prepared tray.

NOTE: Do not interchange measuring spoons and tub lids.

IMPRESSION TAKING

Please take into account the product specific recommendations for the handling of the Aquasil Impression Materials mentioned below. For the light-viscosity materials, please also see the respective directions for use.

APPLICATION	Aquasil Soft Putty	Aquasil Hard Putty	Aquasil Ultra LV/XLV Regular Set	Aquasil Ultra LV/XLV Fast Set
Mixing Time	0:30	0:30	n/a	n/a
Total Time for Application ²	1:30	1:15	2:15	1:15
Minimum Time in Mouth (ISO)	3:00	2:00	2:45	1:45
Setting Time ³	4:30	3:15	5:00	3:00

A) DOUBLE MIX TECHNIQUE

Aquasil Soft Putty or Hard Putty with Aquasil Ultra LV or XLV

Beside common steps for the Double Mix Technique, take into account the special recommendations for the handling of Aquasil Impression Materials mentioned below.

² Describes the time from the start of mix till the impression has to be seated and corresponds to the Total Work Time (ISO) for hand mixed materials and to the Minimum Work Time (including mixing time) (ADA) for cartridge materials.

³ Total time from start of mix to complete set.

Impression Taking

1. Inject light body with cartridge system onto and around tooth preparations. Squeeze dispenser trigger with moderate and even pressure to begin mixing material. Material will stop flowing after trigger is released. It is suggested to syringe completely around preparation at the gingival margins twice, maintaining the syringe tip in the impression material. It is not necessary to blast with air for thinning out the material. Aquasil Ultra LV and XLV Impression Material's hydrophilic properties improve adaptation to the tooth enabling better 'wetting'.
2. Seat tray within 1 1/2 minutes after start of mix of the Regular Set Materials and 1 minute and 15 seconds after start of mix of the Fast Set Materials.
3. Aquasil Impression Materials set faster in the mouth than on the bench. Retain impression in position for at least 3 minutes for the Regular Set Materials and 2 minutes for the Fast Set Materials. Prior to removal, make sure the impression is firm, resilient, and non-tacky.

B) PUTTY/WASH TECHNIQUE

Aquasil Soft Putty or Aquasil Hard Putty with Aquasil Ultra LV or XLV Regular or Fast Set

Due to the excellent flow characteristics it is recommended to choose Aquasil Ultra XLV as the light body for the putty wash technique. Additionally, less material compared to traditional light bodies is used in order to obtain optimal results.

Beside common steps for the Putty/Wash Technique take into account the special recommendations for the handling of Aquasil Impression Materials mentioned below.

Impression Taking

1. For the primary impression seat tray within 1 1/2 minutes after start of mixing Aquasil Soft Putty or 1 minute and 15 seconds after start of mix of Aquasil Hard Putty. Retain impression in position for a minimum of 3 minutes for Aquasil Soft Putty and 2 minutes for Aquasil Hard Putty, respectively.
2. Remove, clean and prepare primary impression for the secondary impression.
3. Clean area to be duplicated with air water spray. Remove spray water.
4. Inject light body with cartridge system onto and around tooth preparations. Squeeze dispenser trigger with moderate and even pressure to begin mixing material. Material will stop flowing after trigger is released. It is suggested to syringe completely around preparation at the gingival margins twice, maintaining the syringe tip in the impression material. It is not necessary to blast with air for thinning out the material. Aquasil Ultra LV and XLV Impression Material's hydrophilic properties improve adaptation to the tooth enabling better 'wetting'.
5. Syringe additional light body on remaining occlusal tooth surfaces and/or on prepared primary impression.
6. Seat tray within 2 minutes and 15 seconds after start of injecting Aquasil Ultra LV/XLV Regular Set and 1 minute and 15 seconds after start of injecting Aquasil Ultra XLV/LV Fast Set, respectively.
7. Aquasil Impression Materials set faster in the mouth than on the bench. Retain impression in position for at least 2 minutes and 45 seconds for Aquasil Ultra LV/XLV Regular Set and 1 minute and 45 seconds for Aquasil Ultra LV/XLV Fast Set. Prior to removal, make sure the impression is firm, resilient, and non-tacky.

CLEANING, DISINFECTION, STORAGE

Rinse impression under cold water and blow dry.

The impression can be disinfected with a standard disinfecting solution or spray.

Do not store impression in direct sunlight.

CASTING AND ELECTROPLATING

To avoid surface porosity of models, delay pouring with plaster for 15 minutes. The impression may be poured up to fourteen days later.

The material is compatible with a range of die materials including epoxy materials and may be silver or copper plated. Wait 60 minutes before start of plating or pouring with epoxy.

TRAY CLEANING INSTRUCTION

Use standard methods to clean the tray after removal of the impression material. Fix Solvent may be used to help removing any remaining layer of Silfix.

STORAGE OF IMPRESSION MATERIAL

Store Aquasil Impression Materials between 10 - 24 °C.

Do not store impression in direct sunlight.

Inadequate storage conditions will shorten the shelf life and may lead to a malfunction of the product.

BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence, which requires identification of the product.

Do not use after expiry date.

RELEVANT MATERIAL PROPERTIES	Aquasil Soft Putty	Aquasil Hard Putty
Delivery system	Tubs	
Mixing technique	Knead with fingertips until colour is uniform	
Mixing ratio	1 : 1	
Classification ISO 4823:1992	Type 0: Very high consistency putty	
Linear dimensional change (ISO)	0.05%	
Maximum strain in compression (ISO)	3.5 %	2.5 %
Recovery from deformation (ISO)	99.5	
Permanent deformation (ADA)	0.5 %	
Detail reproduction	< 50 µm	

If you have any questions, please contact:

Manufacturer:
DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
GERMANY
Phone +49 (0) 75 31 5 83-0

Distributor:
DENTSPLY Limited
Hamm Moor Lane
Addlestone, Weybridge
Surrey KT15 2SE
Phone (0 19 32) 85 34 22

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© DENTSPLY DeTrey 2004-07-06