---

**DENTSPLY**

## Self-Cure Activator

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Magyar

# Self-Cure Activator

**For use in conjunction with XP BOND™ Universal Total-Etch Adhesive or Prime&Bond®NT Nano-Technology Dental Adhesive**

**Self-cure Activator (SCA)** is for use with **XP BOND™** or **Prime&Bond®NT** in conjunction with a dual-cure or self-cure resin cement such as **Calibra®** Esthetic Resin Cement. **SCA** is suitable for bonding of all indirect restorations including metal, porcelain and composite crowns, inlays, onlays, veneers, bridge retainers and endodontic posts. Use of resin cements other than **Calibra** Cement or adhesives other than **XP BOND** or **Prime&Bond NT** with **SCA** is at the discretion and sole responsibility of the dental practitioner.

**XP BOND** and **Prime&Bond NT** are universal self-priming dental adhesives designed for bonding resin-based restorative materials to enamel and dentine.

## COMPOSITION

- Urethane dimethacrylate (UDMA)
- 2-Hydroxyethyl methacrylate (HEMA)
- Catalyst
- Photoinitiators
- Stabilisers
- Acetone
- Water

## INDICATIONS

When used with **XP BOND** or **Prime&Bond NT**, **SCA** is indicated for adhesive cementation of indirect restorations in conjunction with dual-cure or self-cure resin cements.

## CONTRAINDICATIONS

**SCA** and **XP BOND** or **Prime&Bond NT** are contraindicated for:

1. Direct or indirect pulp-capping.
2. Use in patients known to be allergic to methacrylate resins or any of the other components.

## WARNINGS

1. **SCA** contains polymerisable methacrylate monomers. Avoid prolonged or repeated contact with skin (allergic contact dermatitis), oral soft tissues, and eyes. Avoid prolonged inhalation. Do not take internally.
2. **Eye Contact:** **SCA** contains methacrylates which may be irritating to eyes. Before using this product wear protective glasses and cover the patient's eyes to protect from splashing material. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
3. **Skin Contact:** **SCA** contains polymerisable monomers which can cause skin sensitisation (allergic contact dermatitis) in susceptible individuals. In case of contact with skin, immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitisation or other allergic reactions occur, discontinue use and seek medical attention.
4. **Oral mucosa contact:** Avoid contact with oral soft tissues. If accidental contact occurs, flush mucosa with plenty of water and expectorate water. If sensitisation of mucosa persists, seek medical attention immediately.

5. Ingestion: Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. If nausea or illness develops, seek medical attention immediately.
6. SCA contains acetone. Do not inhale vapor. Inhalation: Give oxygen or artificial respiration if necessary.

## PRECAUTIONS

1. This product is only intended for use as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the dental practitioner.
2. SCA is for use in conjunction with XP *BOND* or Prime&Bond *NT* and DeTREY® Conditioner 36 For correct use, please also refer to the Directions for Use of the respective products.
3. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
4. Dental materials containing eugenol should not be used in conjunction with this product because they may impede curing and cause softening of the polymeric components of the material.
5. Insufficient data exist to support the use of desensitising agents and/or cavity cleansing agents with SCA. Recommended pre-treatment is pumice or prophylaxis paste with a rubber cup.
6. Contact with saliva, blood and/or some astringent solutions during adhesive procedures may cause failure of the restoration. Use of a rubber dam or adequate isolation is recommended.
7. Exercise caution if mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures. In-vitro studies have suggested compromised adhesive bond strength to enamel and dentine contaminated with haemostatic agents. Marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.
8. In-vitro studies have shown a relationship between adhesive performance and airthinning/drying/solvent evaporation techniques. Improper over-thinning or under-evaporating may compromise bond strength, leading to microleakage, post-operative sensitivity and/or restoration failure. (See STEP-BY-STEP INSTRUCTIONS)
9. The use of a dual-cure adhesive system, such as XP *BOND* or Prime&Bond *NT*, with SCA can shorten the working time of a dual-cure resin cement system. This effect should be investigated in the laboratory prior to clinical use.
10. The SCA bottle should be tightly closed immediately after use.
11. When mixed with SCA, XP *BOND* and Prime&Bond *NT* can still be cured with visible light. Proceed immediately once materials have been placed in mixing well or protect from ambient light.
12. Must only be used in well ventilated areas.
13. Flammable: SCA contains acetone. Keep away from sources of ignition.
14. **Storage:** SCA should be kept out of direct sunlight and stored in a well ventilated place at room temperature not exceeding 25 °C. Refrigerated storage is not required, but is acceptable when not in use. Inadequate storage conditions will shorten the shelf life and may lead to the product malfunctioning. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

## ADVERSE REACTIONS

1. Product may irritate the eyes and skin.  
Eye contact: Irritation and possible corneal damage.  
Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous membranes: Inflammation, oedema, sloughing. (See WARNINGS)
2. Product may affect health seriously if ingested. (See WARNINGS)

3. Inhalation of vapors may cause varying degrees of damage to the affected tissue and also increased susceptibility to respiratory illness. (See PRECAUTIONS)
4. The following medical conditions are generally aggravated by exposure: Individuals with preexisting skin disorders, eye problems or impaired liver, kidney and respiratory functions.

## STEP-BY-STEP INSTRUCTIONS

### General Information

In conjunction with Prime&Bond *NT* or *XP BOND*, SCA is indicated for cases where application of a dual-cure or self-cure resin cement is intended. SCA will promote adhesion of compatible self-cure or dual-cure resin cements to the adhesive layer.

In cases where it would be impossible or undesirable to light-cure the adhesive layer, SCA is to be used in conjunction with *XP BOND* and mixed Calibra (base + catalyst). Additional lightcuring of the *XP BOND*/SCA adhesive layer is optional.

When using SCA with Prime&Bond *NT*, light-curing the adhesive layer is mandatory.

To ensure that indirect restorations/post surfaces are pretreated correctly (e.g. microroughening, silanization with Calibra Silane Coupling Agent), please refer to the Directions for Use of the respective products. To ensure fast luting, it is advisable to prepare the restoration prior to applying the adhesive.

### Applying SCA with Prime&Bond *NT* or *XP BOND*

1. Pretreat the prepared tooth/cavity by cleaning, protecting the pulp and acid-etching as described in the Prime&Bond *NT* or *XP BOND* Directions for Use.
2. Dispense Prime&Bond *NT* or *XP BOND* first, followed by SCA, at a ratio of 1:1 into a suitable dappen dish (e.g. Apli-Dish™).
3. Mix liquids for at least 1 or 2 seconds with a standard applicator tip.
4. If (as in the case with Calibra) mixed adhesive is applied onto an internal bonding surface of the restoration, air dry for at least 5 seconds and light-cure for 10 seconds<sup>1</sup>. Light-curing is optional when *XP BOND*/SCA are used in conjunction with mixed Calibra (base + catalyst).
5. Immediately apply adhesive mixture to the cavity. Wet all cavity surfaces uniformly. With Prime&Bond *NT*/SCA, apply ample amounts of mixture to the cavity. With *XP BOND*/SCA, avoid pooling.
6. Leave the surfaces undisturbed for 20 seconds.
7. Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds. Surface should have a uniform, glossy appearance. If not, repeat steps 5 and 7. Surface should not show areas of excessive adhesive thickness or pooling. Repeat air drying/evaporation outlined above if necessary.
8. Light-cure adhesive layer for a minimum of 10 seconds<sup>1</sup>. Ensure uniform exposure of all cavity surfaces. Light-curing is optional when *XP BOND*/SCA are used in conjunction with mixed Calibra (base + catalyst).

Immediately apply compatible resin-based luting cement in desired curing mode, i.e. self-cure (e.g. for PFM crowns, metal posts) or dual-cure (e.g. for indirect porcelain/composite inlays, translucent fibre posts). Seat indirect restoration immediately. For further information on how to use Calibra resin cement or other resin based luting cements, please refer to the Directions for Use of the respective product.

## CLEANING AND DISINFECTION

To avoid splattering or spraying the Self-cure Activator bottle with body fluids or placing it in contact with contaminated hands or oral tissues, use of a protective barrier is recommended. Repeated disinfection may damage label.

## BATCH NUMBER AND EXPIRY DATE

Do not use after expiry date. ISO standard is used: "YYYY/MM"

The batch number, the reorder number and the expiry date should be quoted in all correspondence which requires identification of the product.

If you have any questions, please contact:

### EU-Authorized Representative:

**DENTSPLY DeTrey GmbH**

**De-Trey-Str. 1**

**78467 Konstanz**

**GERMANY**

**Phone +49 (0) 75 31 5 83-0**

### Distributor:

**DENTSPLY Limited**

**Hamm Moor Lane**

**Addlestone, Weybridge**

**Surrey KT15 2SE**

**Phone (0 19 32) 85 34 22**

©2006-04-03 DENTSPLY International Inc. All Rights Reserved.

---

<sup>1</sup> Check your curing light for a minimum output level of 800 mW/cm<sup>2</sup>. Light cure adhesive for at least 20 seconds if the light intensity is between 550 and 800 mW/cm<sup>2</sup>.

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