DENTSPLY

Aquasil

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Česky

Aquasil *Bite*

Vinyl Polysiloxane Bite Registration Material Cartridge Delivery System

CAUTION: For dental use only. In the U.S., federal law restricts this device to sale by or on the order of a dentist.

DESCRIPTION

Aquasil Bite is designed for accurate interocclusal and jaw relation records. Aquasil Bite offers minimum resistance to closure during the 30 seconds working time, minimizing deviation of the mandible producing consistently accurate jaw relation records. Aquasil Bite accurately records dental morphology, enabling the dentist and lab technician to easily articulate models without interferences from convex tooth surfaces. The rapid rigidity developed at 1 minute mouth removal time reduces inaccuracies created by mandibular movement during setting and allows accurate model articulation with a minimal amount of distortion producing elastic rebound.

The accurate detail reproduction, minimal distortion and quick set also allow use of Aquasil Bite as a pre-preparation quadrant matrix material in the direct fabrication technique of provisional (temporary) restorations.

Aquasil Bite is available in 50 ml cartridge delivery.

COMPOSITION

Aquasil Bite is an addition-reaction silicone (vinyl polysiloxane) composed of:

- Polydimethyldivinylsiloxane polymer
- Silicon Dioxide
- Polymethylhydrogen Siloxane
- Surfactant
- Plasticizer
- Pigments

INDICATIONS FOR USE

Intraoral recording of jaw relations and occlusal relationships, used in laboratory cast mounting procedures. May be used to make registrations between natural or restored teeth and/or between complete or partially edentulous occlusion rims.

Pre-operative quadrant matrix impression for use in the direct fabrication of provisional (temporary) crowns, short-span bridges, inlays and onlays.

CONTRAINDICATIONS

None known.

TECHNICAL DATA

Work time	30 seconds minimum from the start of mix at 22 °C
Minimum Removal Time (MRT)	1 minute from placement onto occlusal surfaces
% Recovery from deformation	> 99 %
% Strain in compression	< 1 %
Detail Reproduction	2 μ
Linear dimension change	0.05 %
Mixing ratio by volume	1 part base to 1 part catalyst
Shore-A hardness	95

WARNINGS

1. Avoid prolonged or repeated exposure of the Aquasil Bite with skin and eyes. Irritation and possible corneal damage may result. Skin rash, oral mucosa irritation, or other allergic reactions may result in susceptible individuals.

Eye and Skin Contact: Flush eyes with flowing water for 15 minutes and consult a physician. Flush skin with flowing water for 15 minutes, then wash area with soap and water after contact. Flush oral tissues with copious amounts of water. Consult a physician if rash persists.

Ingestion: Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. This material is not hazardous when small quantities are ingested. Larger quantities may cause bowel obstruction. Seek medical attention in the event of digestive irregularities.

2. Do not use Aquasil Bite as a temporary reliner. Use of the material in this manner may cause irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.

PRECAUTIONS

1. This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
2. Allow Aquasil Bite to reach room temperature prior to use. Higher temperatures reduce work times and laboratory bench set times (faster), lower temperatures increase them (slower). Clinically, intraoral set time (Mouth Removal Time) is unaffected within the range of allowable storage and initial mixing temperatures.
3. Due to the stiffness of Aquasil Bite, severe undercuts should be blocked out prior to use.
4. Material should extrude easily. DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause cartridge rupture.
5. Insufficient data exist to support the use of Aquasil Bite as a final impression material for definitive cast, appliance, restoration or prosthesis fabrication.
6. Insufficient data exist to support the use of Aquasil Bite for registrations on exposed bone. Aquasil Bite should only be used on or around intact mucosa.

ADVERSE REACTIONS

1. Corneal damage may result with prolonged eye exposure to the registration material. See Warnings.
2. Allergic contact dermatitis, oral mucosa irritation, and other allergic reactions may occur in susceptible individuals. See Warnings.
3. Skin irritation may result from direct exposure and skin defatting from prolonged exposure. See Warnings and Precautions.
4. Bowel obstruction or other digestive distress may result from ingestion of mixed registration material. See Warnings.

INTERACTIONS

1. Do not contaminate Aquasil Bite or teeth/surfaces to be recorded with latex or "sulfur-containing" polymeric gloves. Residue from some brands of latex gloves may leach out especially in the presence of blood, saliva and some astringents and retraction cords. These residues on the surfaces of teeth can interfere with the chemical set of the polyvinyl surface and lead to loss of surface detail and inaccuracies. Non-latex gloves worn during manipulation will avoid this possibility.
2. Use of some hand soaps and lotions can lead to interference with setting reaction. Contamination may accumulate on the surface of the skin and may not easily wash off. Non-latex gloves worn during manipulation will avoid this possibility.

STEP-BY-STEP INSTRUCTIONS

Cartridge Dispensing Gun Instructions

1. Raise the release lever vertically upward, while simultaneously pulling the plunger all the way back in the dispenser handle.
 2. Cartridge Loading
 - 2.1 Open cartridge lock by lifting up top clasp.
 - 2.2 Orient and insert cartridge with v-shaped notch facing down.
 - 2.3 Close the top clasp to lock cartridge into dispenser gun.
 - 2.4 Remove the cartridge cap by turning 90° counterclockwise. Discard the cartridge cap. The used mixing tip can be left in place until the next use after disinfection to serve as a self-sealing cap.
 - 2.5 Dispense a small amount of base and catalyst before installing mix tip to ensure even flow from cartridge. Use gentle pressure. Be sure no plug is present. If any plug blocks flow, clear with an instrument. Wipe away excess from the cartridge carefully so base and catalyst do not cross contaminate and cause obstruction of the nozzle.
 - 2.6 Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside colored rim of the mix tip with the v-shape notch on the cartridge flange. If mix tip does not seat easily, be sure the internal end view of the mix tip is aligned correctly. The mix tip itself can turn in its cap. The two holes should line up with the v-shaped piece that is on the outside rim of the mix tip cap to facilitate inserting. Turn tip to align in its cap.
 - 2.7 When mix tip colored rim v-notch is lined up with the v-shaped notch on the cartridge flange, turn colored mix tip cap 90° in a clockwise direction to lock in place on cartridge.
 3. Dispenser Operation
 - 3.1 Squeeze dispenser trigger with moderate and even pressure to begin mixing material. Material will stop flowing after trigger is released.
 - 3.2 The Aquasil Bite may be dispensed directly to the mouth with provided mix tips.
- OR**
- An intraoral syringe may be filled by extruding the mixed material directly into the syringe by back filling.
- 3.3 Disinfect cartridge dispensing gun, cartridge and attached used mix tip for storage.

Immediately prior to subsequent usage, remove the used mixing tip and squeeze the dispensing gun trigger to dispense (bleed) material from the two openings of cartridge. Wipe end of cartridge. Install a new mixing tip and twist one-quarter turn to lock to position. To remove the cartridge, hold the release lever vertically while simultaneously withdrawing plunger in a backward direction. Dispose of cartridge when empty. Do not attempt to refill or consolidate used cartridges.

1. Interocclusal Bite Registration

Dispensing Material

1. Inject material directly onto the patient's occlusal surfaces and/or occlusion rims. The areas where teeth have been prepared should be injected first. Usually, only teeth of the upper or lower arch need be covered. However, in recording the anterior teeth, it is suggested that material be placed over the teeth of both arches. A layer of about 5 mm of Aquasil Bite is usually sufficient. Do not exceed 30" while injecting.

Technique Tip: This material may be used with or without a bite registration tray, but in most applications, a tray will not be required. In those situations where use of a tray is desired, the mixed material is placed on both sides of the tray to a depth of 2-3 mm. This material has been found to work well with mesh fabric type trays.

2. Gently guide or instruct patient to close into desired position. Due to the thixotropic nature, Aquasil Bite offers no resistance to closure during the 30 seconds work time.
3. Instruct the patient to maintain, or manually maintain jaw relationship without movement until 1 minute from end of placement. Check material in the mouth (not on the bench) to be sure it is firm, resilient, and non-tacky before removal.

Removing Aquasil Bite

1. After 1 minute from end of placement, the material may be removed from the patient's mouth. Have the patient open and carefully release the bite registration from the teeth. Check for residual material in undercut areas.
2. Rinse registration under cold water and blow dry.

2. Matrix for fabrication of provisional (temporary) restorations

Consult manufacturer's complete Directions for Use of selected provisional material for vinyl polysiloxane matrix technique compatibility.

Prior to preparation of the tooth/teeth, evaluate existing contours and anatomy. To duplicate existing anatomy proceed to matrix impression. If contours or anatomy are deficient, light cured composite may be placed and cured without bonding agent on the tooth/teeth as a mock up to approximate the desired contours. A missing tooth may be temporarily replaced with an acrylic denture tooth prior to matrix impression.

Matrix Impression

1. Select and/or prepare a suitable quadrant metal tray, a firm quadrant disposable tray or a quadrant custom tray. Use rigid trays of sufficient size to provide at least 2-3 mm thickness of matrix impression material.
2. Dispense Aquasil Bite Registration Material directly into impression tray from cartridge. Do not layer the material or place it into rows. This technique ensures the least amount of air incorporation.

Technique Tip: Loaded tray must be inserted within 30" (thirty seconds) from time of first dispensing.

3. Seat loaded tray. Tray must be seated within 30" (thirty seconds). Retain matrix impression in position until firmly set. (**NOTE:** The material sets faster in the mouth than on the bench). Check matrix impression in the mouth (not on the bench) to be sure it is firm, resilient, and non-tacky before removal. The Minimum Removal Time is 1 minute from the start of mix. At 22 °C Aquasil Bite Registration Material has a minimum work time of 30 seconds. Variations in temperature from that stated above will affect work and set times. Higher temperatures reduce work times and lower temperatures increase them.
4. Remove matrix impression by pulling slowly to break seal, then snap out along the long axis of the teeth. Check for residual material in undercut areas. Rinse matrix impression under cold water and blow dry.

5. Complete preparation(s) and fabricate provisional (temporary) restorations as per provisional material manufacturer's Directions for Use.

CLEANING

Dispenser and cartridge cleaning instructions.

The cartridge dispensing guns may be cleaned by scrubbing with hot water and soap or detergent.

DISINFECTION

Dispenser and cartridge disinfection instructions

Dispensers, cartridges and used tips left in place on cartridges for storage, exposed to splatter or spray of body fluids or that may have been touched by contaminated hands, or oral tissues, should be disinfected with a hospital-level disinfectant. Acceptable disinfectants are those that are registered as tuberculocidal. Iodophors, sodium hypochlorite (5.25 %), chlorine dioxide and dual or synergized quaternary ammoniums are approved disinfectants. Water-based disinfectant solutions are preferred. Use impregnated wipes not sprays.

Note: Wipe cartridge GENTLY, vigorous wiping may destroy the label.

Some phenolic based agents and iodophor based products may cause surface staining.

Agents containing organic solvents, such as alcohol may tend to dissolve the plastic dispenser and plunger. The disinfectant manufacturer's directions should be followed properly for optimum results.

Technique Tip: Unset Aquasil Bite can be removed from clothing with dry cleaning solvent.

Disinfection of the Bite Registration Material

The registration should be disinfected with a hospital-level disinfectant. Acceptable disinfectants are those that are registered as tuberculocidal. Iodophors, sodium hypochlorite (5.25 %), chlorine dioxide, and dual or synergized quaternary ammoniums are approved disinfectants. To disinfect polyvinyl siloxane material, thoroughly soak by spraying or immerse the registration in any recommended hospital level disinfectant except neutral glutaraldehyde for the contact time recommended by the disinfectant manufacturer for optimum results. Water-based disinfectant solutions are preferred.

MODEL PLACEMENT INTO REGISTRATION

The registration should be removed from the disinfectant, rinsed with water and exposed to air to dry for at least the same amount of time it was exposed to the disinfectant solution before it is used. Do not store registration in direct sunlight. If the registration is to be shipped, use suitable packaging to preclude distortion.

Any excess material extending onto soft tissue areas, into interproximal areas or gingival to the height of contour of the teeth should be removed to obtain maximum accuracy when seating the registration on the model. Best results are obtained by trimming with a sharp scalpel blade or fine, sharp scissors.

STORAGE

Store Aquasil Bite between 2 and 28 °C. Protect from moisture. Do not freeze. Do not store used cartridges with new (unused) mix tips attached.

Inadequate storage conditions may shorten the shelf life and lead to malfunction of the product.

BATCH NUMBER () AND EXPIRY DATE ()

1. Do not use after expiration date. ISO standard uses: "YYYY/MM."
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number on the cartridge
 - Expiration date

If you have any questions, please contact:

Manufacturer ():
DENTSPLY DeTREY GmbH
De-Trey-Str. 1
78467 Konstanz
Germany
Phone +49 (0) 75 31 5 83-0

Distributor:
DENTSPLY Limited
Hamm Moor Lane
Addlestone, Weybridge
Surrey KT15 2SE
Phone (0 19 32) 85 34 22

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© DENTSPLY DeTREY 2007-02-01