# Aquasil Putty DECA™

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# Aquasil Putty DECA™

## Vinyl Polysiloxane Impression Material for Dynamic Mixing

CAUTION: For dental use only.  
USA: Rx only

Content	Page
<b>1 Product Description</b> .....	<b>2</b>
<b>2 General Safety Notes</b> .....	<b>3</b>
<b>3 Step-by-step Instructions</b> .....	<b>5</b>
<b>4 Hygiene</b> .....	<b>7</b>
<b>5 Lot number and expiration date</b> .....	<b>8</b>

## 1 Product Description

Aquasil Putty DECA™ is a **hydrophilic base impression material** consisting of addition curing vinyl silicones. The mixing ratio is 5 volumes base paste and 1 volume catalyst paste.

Aquasil Putty DECA™ cartridges can be used with suitable mixing machines.

Aquasil Putty DECA™ shows **hydrophilic properties** and is compatible with the hydrophilic Aquasil® Ultra syringe materials and moist surfaces. Therefore it yields an exact reproduction of subgingival areas and other moist surfaces.

Aquasil Putty DECA™ meets the requirements of **ISO 4823** (2000), **type 0** (putty consistency).

### 1.1 Delivery forms

Aquasil Putty DECA™ is available in 5:1 DECA™ delivery.

### 1.2 Composition

- Vinyl siloxane
- Hydrogen siloxane
- Silica
- Paraffine

### 1.3 Indications

- Aquasil Putty DECA™ is designed as a firm tray material suitable for: all tooth, implant (closed flap) or tissue impression techniques.
- Aquasil Putty DECA™ is a very high viscosity type impression material with ideal mechanical properties as tray material for the putty/wash technique.
- It may also be used as tray material for the Dual phase one step quadrant arch, double arch or full arch impression technique (double mix technique) with Aquasil® Ultra LV or XLV Impression Materials.

Classification ISO 4823	Aquasil Putty DECA™: Type 0: Putty consistency
<b>Total Working Time (23 °C/74 °F)</b>	<b>1'30"</b>
<b>Recommended time in mouth</b>	<b>3'00"</b>
Polymerisation shrinkage after 24 h	0.05-0.3%
Elastic Recovery	≥ 98.8%
Strain in compression	2.3-3.3%
Mixing ratio by volume (DECA™ delivery)	5 parts base to 1 part catalyst

Data on File.

## 1.4 Contraindications

None known.

## 2 General Safety Notes

Be aware of the following general safety notes and the special safety notes in other chapters of these instructions.

### 2.1 Warnings

1. Avoid prolonged or repeated exposure of Aquasil Putty DECA™ with skin and eyes. Irritation and possible corneal damage may result. Skin rash, oral mucosa irritation, or other allergic reactions may result in susceptible individuals. Do not take internally.  
 Eye and Skin Contact: Flush eyes with flowing water for 15 minutes and consult a physician. Flush skin with flowing water for 15 minutes, then wash area with soap and water after contact. Flush oral tissues with copious amounts of water. Consult a physician if rash persists.  
 Ingestion: Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. This material is not hazardous when small quantities are ingested. Larger quantities may cause bowel obstruction. Seek medical attention in the event of digestive irregularities.
2. Do not use Aquasil Putty DECA™ as a temporary reliner. Use of the material in this manner may cause irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.

### 2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.
2. Aquasil Putty DECA™ should not be intermixed with, or used in conjunction with, other manufacturer's vinyl polysiloxanes. Aquasil Putty DECA™ should not be intermixed with, or used in conjunction with, polyether, polysulfide or conventional (condensation cured) silicones.
3. Insufficient data exist to support the use of Aquasil Putty DECA™ for impressions on exposed bone. Aquasil Putty DECA™ should only be used on or around intact mucosa.
4. Allow Aquasil Putty DECA™ to reach room temperature prior to use. Higher temperatures reduce work times and laboratory bench set times (faster), lower temperatures increase them (slower). Clinically, intraoral set time (Mouth Removal Time) is unaffected within the range of allowable storage and initial mixing temperatures.

5. To obtain optimum physical properties, the mix should be completely homogeneous (streak free).
6. Do not contaminate retraction cords with latex or “sulfur-containing” polymeric gloves. Handle cords with college pliers. Do not touch prepared tooth with gloves.
7. Use of some hand soaps and lotion can lead to interference with setting reaction. Contamination may accumulate on the surface of the skin and may not easily wash off. Non-latex gloves worn during mixing will avoid this possibility.
8. If the preparation is treated with hydrogen peroxide, rinse with copious amounts of water before making the impression.
9. Aquasil Putty DECA™ should be used with suitable impression trays pretreated with tray adhesive (**see Step-by-Step Instructions**). Refer to tray adhesive manufacturer’s Instructions for Use for compatibility.
10. Material should extrude easily: **DO NOT USE EXCESSIVE FORCE**. Excessive pressure may result in unanticipated extrusion of the material or cause cylinder rupture.

### 2.3 Storage conditions

Store Aquasil Putty DECA™ at temperatures between 10-25 °C/50-77 °F. Allow material to reach room temperature prior to use. Working time is affected by temperature (**see Precautions, Step-by-Step Instructions**). Protect from moisture. Do not freeze. Do not use after expiration date.

### 2.4 Adverse Reactions

1. Corneal damage may result with prolonged eye exposure to the impression material (**see Warnings**).
2. Allergic contact dermatitis, oral mucosa irritation, and other allergic reactions may occur in susceptible individuals (**see Warnings and Precautions**).
3. Skin irritation may result from direct exposure and skin defatting from prolonged exposure (**see Warnings and Precautions**).
4. Bowel obstruction or other digestive distress may result from ingestion of mixed impression material (**see Warnings**).

### 3 Step-by-Step Instructions

#### Aquasil Putty DECA™ Loading Instructions

The cylinder fits the Duomix<sup>1</sup> and Pentamix 2<sup>1</sup> dynamic mixers. Carefully remove dispensing cylinders from the package.

**Note:** The cartridges cylinders are heavy and dropping may result in damage to the cylinders rendering it unusable.

**DECA™ Cylinder Loading:** Also refer to separate dynamic mixer manufacturer's instructions.

1. Firmly hold and support the full cylinder with dispensing nozzles up.
2. Grasp the tab on the cap in the dispensing nozzles and bend in an upward direction 90°. Do not place undue stress on tab.
3. Firmly grasp the tab with thumb and forefinger while supporting the cylinder and remove the cap from the cylinder. There will be some resistance to cap removal.  
**Note:** Do not drop the cylinder. The cylinder may break rendering it unusable.
4. Insert the cylinder into the dynamic mixer according to the instructions supplied by the manufacturer.
5. Following the machine manufacturer's instructions to dispense, press control button. Bleed the material from the two nozzles of cartridge. (This should be done prior to each use.) Wipe off the extruded material in a vertical direction to avoid cross contamination of the nozzles.  
**Note:** Discard the first 3 cm/1.18 in of impression material extruded from newly filled cartridges.
6. Install a supplied mixing tip on the extrusion nozzles following the machine manufacturer's instructions. (**Note:** Do not use other manufacturer's mix tips – they will not fit snugly and can cause material to leak from extrusion nozzles.)
7. If the mixing tip does not seat easily, be sure the internal central hex of the dynamic mix tip is aligned correctly with the drive shaft of the dynamic mixer.
8. Once the mixing tip is properly placed, place yellow nozzle locking ring over the mix tip and turn to the right (1/4 turn) until ring is locked in place.
9. Check whether the hexagonal drift shaft does engage when starting. The spiral in the mixing tip must rotate when starting.
10. Dispense the mixed base and catalyst into a suitable tray.
11. Immediately after dispensing, back the plungers of the machine off the cylinder to prevent continued pressure on the material and extrusion of excess material from the tip.
12. After disinfection, keep the used mixing tip in place until the next use to serve as a self-sealing cap. To remove the used mixing tip turn the bayonet ring anti-clockwise and draw it off.

#### Preparation of syringe/wash material

If using dual-phase technique, suitable Aquasil® Ultra syringe (wash) impression material in the cartridge or digit® Targeted Delivery System Unit Dose Dispensed delivery (available separately) may be selected and prepared for use according to the Instructions for Use.

---

<sup>1</sup> Not registered trademarks of Dentsply Sirona, Inc.

## Impression Techniques

### 3.1 Putty/wash technique

#### **Aquasil Putty DECA™ with Aquasil® Ultra LV or XLV Regular or Fast Set**

Due to the excellent flow characteristics it is recommended to choose Aquasil® Ultra XLV as the light body for the putty wash technique. Additionally, less material compared to traditional light bodies is used in order to obtain optimal results.

Beside common steps for the putty/wash technique take into account the special recommendations for the handling of Aquasil® Impression Materials mentioned below.

#### **Impression Taking**

1. Dispense Aquasil Putty DECA™ directly into impression tray from dynamic mixer. For the primary impression seat tray not later than **1 and 1/2 minutes** (working time at 23 °C or 74 °F) after extruding Aquasil Putty DECA™ onto the tray. Retain impression in position for a minimum of 3 minutes.
2. Remove, clean and prepare primary impression for the secondary impression.
3. Clean area to be duplicated with air water spray. Remove spray water.
4. Inject light body with cartridge system onto and around tooth preparations. Squeeze dispenser trigger with moderate and even pressure to begin mixing material. Material will stop flowing after trigger is released. It is suggested to syringe completely around preparation at the gingival margins twice, maintaining the syringe tip in the impression material. It is not necessary to blast with air for thinning out the material. Aquasil® Ultra LV and XLV Impression Material's hydrophilic properties improve adaptation to the tooth enabling better 'wetting'.
5. Syringe additional light body on remaining occlusal tooth surfaces and/or on prepared primary impression.
6. Seat tray within 2 minutes and 15 seconds after start of injecting Aquasil® Ultra LV/XLV Regular Set and 1 minute and 15 seconds after start of injecting Aquasil® Ultra LV/XLV Fast Set, respectively.
7. Aquasil® Impression Materials set faster in the mouth than on the bench. Retain impression in position for 5 minutes from start of mix for Aquasil® Ultra LV/XLV Regular Set and 3 minutes from start of mix for Aquasil® Ultra LV/XLV Fast Set. Prior to removal, make sure the impression is firm, resilient, and non-tacky.

The impression should be stored under normal room conditions (max. 25 °C or 78 °F). Used mixing tip will be left in place until next use.

### 3.2 Dual Phase one-step Technique: with a Rigid Tray or with Double Arch Technique.

**Technique Tip:** When using a dynamic mixer, it is highly recommended to utilize an assistant to assure simultaneous syringing and tray loading.

1. Select and/or prepare a suitable metal tray, a firm disposable tray or a custom tray. Use rigid trays of sufficient size to provide at least 2-3 mm thickness of impression material.
2. Brush a thin layer of tray adhesive following the tray adhesive manufacturer's Instructions for Use. Adhesive should be used even with perforated trays.
3. If there is any debris or blood, clean the field with air water spray.
4. Remove excess water with suction. Do not desiccate the tooth. Dry and isolate with cotton rolls.

#### **For syringe material:**

5. Attach an intraoral tip to the end of the mix tip for direct intraoral syringing from a cartridge or prepare digit® Targeted Delivery System Unit Dose Dispensed syringe. To simultaneously inject the syringe material and load the tray will require a chairside assistant.

#### **For tray material:**

- 6a. Dispense Aquasil Putty DECA™ material (assistant) directly into impression tray from the dynamic mixer. Do not layer the material or place it into rows. This technique ensures the least amount of air incorporation.
- 6b. Remove cotton rolls. It is suggested to simultaneously syringe (dentist) completely around the margins of the preparation/abutment, making sure to keep intraoral tip within impression material. Inject syringe material into any existing anatomy and continue syringing around the preparation(s)/abutment(s) until the preparation(s)/abutment(s) is completely covered with syringe material. Do not blast with air. Aquasil® Ultra Impression Material's hydrophilic properties improve adaptation to the tooth enabling better "wetting".
7. Seat loaded tray within 1 minutes and 10 seconds from time of first introduction of Aquasil® Ultra LV/XLV Regular Set intraorally and 35 seconds from time of first introduction of Aquasil® Ultra LV/XLV Fast Set intraorally. Retain impression in position until firmly set. (**Note:** The material sets faster in the mouth than on the bench.) Check impression in the mouth (not on the bench) to be sure it is firm, resilient, and non-tacky before removal. The Minimum Removal Time is 5 minutes from the start of mix when used with Aquasil® Ultra LV/XLV Regular Set, 4 minutes and 30 seconds from the start of mix when used with Aquasil® Ultra LV/XLV Fast Set. Higher temperatures reduce work times and lower temperatures increase them.
8. Remove impression by pulling slowly to break seal, then snap out along the long axis of the tooth. Rinse impression under cold water and blow dry.

## **4 Hygiene**

### **4.1 Cleaning**

The cylinders may be cleaned by scrubbing with hot water and soap or detergent. Excess material may be cleaned with alcohol-moistened gauze.

### **4.2 Disinfection**

Cylinders and used tips left in place on cylinders for storage should be disinfected with a hospital-level tuberculocidal disinfectant according to national and local regulations. Use impregnated wipes not sprays.

**Note:** Wipe cylinder GENTLY, vigorous wiping may destroy the label. Some phenolic-based agents and iodophor-based products may cause surface staining. Agents containing organic solvents, such as alcohol may tend to dissolve the plastic dispenser and plunger. The disinfectant manufacturer's instructions should be followed properly for optimum results.

**Technique Tip:** Unset Aquasil® Impression Material can be removed from clothing with dry cleaning solvent.

Clean and disinfect Dynamic Mixing Machine following machine manufacturer's instructions.

### **4.3 Disinfection of the Impression Material**

The impression should be disinfected with a hospital-level tuberculocidal disinfectant according to national and local regulations. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide, and dual or synergized quaternary ammoniums are approved disinfectants. To disinfect polyvinyl siloxane material, thoroughly soak by spraying or immerse the impression in any recommended hospital level disinfectant except neutral glutaraldehyde for the contact time recommended by the disinfectant manufacturer for optimum results. Water-based disinfectant solutions are preferred.

#### 4.4 Casting and Electroplating

The impression should be removed from the disinfectant, rinsed with water and exposed to air to dry for at least the same amount of time it was exposed to the disinfectant solution before it is poured. The impression can be poured 30 minutes after removal from the mouth and up to fourteen days.

Impressions may be silver- or copperplated.

Do not store impression in direct sunlight. If the impression is to be shipped, use suitable packaging to preclude distortion. The material is compatible with a range of die stones.

#### 5 Lot number ( ) and expiration date ( )

1. Do not use after expiration date.  
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
  - Reorder number
  - Lot number
  - Expiration date

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© Dentsply France 2008-06-26
Dentsply France SAS