Dycal[®]

Radiopaque Calcium Hydroxide Composition

| | |
|--------------------------------|----|
| Directions for Use _____ | 3 |
| Gebrauchsanweisung _____ | 9 |
| Mode d'emploi _____ | 15 |
| Gebruiksaanwijzing _____ | 21 |
| Istruzioni per l'uso _____ | 27 |
| Instrucciones de uso _____ | 33 |
| Instruções de utilização _____ | 39 |
| Οδηγίες χρήσης _____ | 45 |
| Bruksanvisning _____ | 51 |
| Brugsanvisning _____ | 57 |
| Инструкция по применению _____ | 63 |
| Instrukcja użytkowania _____ | 69 |
| Návod k použití _____ | 75 |

Dycal[®]

Radiopaque Calcium Hydroxide Composition

Directions for Use - English

Dycal® is a rigid, self-setting, radio-opaque calcium hydroxide composition. It will not inhibit the polymerization of acrylic and composite restorations. Dycal is available in tubes and in CFC-free pressure packs.

COMPOSITION

Base paste

1,3-Butylene glycol disalicylate
Zinc oxide
Calcium phosphate
Calcium tungstate
Iron oxide pigments

Catalyst paste

Calcium hydroxide
N-ethyl-o/p-toluene sulfonamide
Zinc oxide
Titanium dioxide
Zinc stearate
Iron oxide pigments (dentine shade only)

INDICATIONS

Direct and indirect pulp capping.
Protective liner under dental filling materials, cements and other base materials.

CONTRAINDICATIONS

Hypersensitivity against sulfonamides or other components, respectively, of the product.

SIDE EFFECTS

May cause allergic reactions in persons susceptible to sulfonamides.

NEGATIVE INTERACTIONS WITH OTHER DENTAL PRODUCTS

Acetone-containing products are detrimental to calcium hydroxide materials. Minimize contact of these products with calcium hydroxide liners.

APPLICATION

Dosage

Extrude equal volumes of the base and catalyst paste (1.17 g : 1.00 g) on the mixing pad provided. Using the pressure pack, this is achieved by applying light sideways pressure to the tip on the mixing pad of the container

nozzle. The nozzle is empty before the first use and it may therefore take a few seconds before the material is extruded.

Do not attempt to control setting time by increasing or decreasing the amount of catalyst dispensed. The product is adjusted to provide a constant working time in spite of such variations.

Mixing

Stir immediately using a Dycal® applicator until a uniform colour is achieved. Do not overspatulate. To provide maximum working time, complete mixing of the two pastes within 10 seconds.

Application

A ball-pointed instrument like the Dycal® applicator is an

effective instrument for placing the mix in the cavity and spreading it over the floor to the desired depth.

Quickly apply the mixed material into the cavity to take advantage of the essential flowability of the mass for efficient placement before setting starts. The cavity preparation should be dry.

Avoid touching the margins of the cavity.

Avoid placing a large bulk of material.

Setting

The mixed material will set in 2½ to 3½ minutes on the mixing pad (23 °C with 50% relative humidity). The set time is shorter in the mouth due to moisture and temperature. Increased humidity and elevated

temperature will reduce available work time and set time. Conversely, decreased humidity and lower temperature will increase available work time and set time.

Remove any excess set material from retention areas or margins with a sharp spoon excavator or bur.

STORAGE

Keep dispensers or tubes tightly closed.

Store at room temperature or below. Do not pierce or burn pressure packs even when empty.

Separation of solid and liquid phases of the base paste is accelerated at elevated temperatures. When stored at reduced temperatures, allow the material to return to normal room temperature before using.

BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

SPECIAL NOTES

1. For dental use only.
2. Dycal catalyst paste contains calcium hydroxide. Avoid contact with skin and eyes. If accidental contact occurs, flush area with generous amounts of water. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
3. Dycal base paste contains 1,3-butylene glycol disalicylate which is irritating to skin and eyes. In case of contact with eyes, rinse immediately with plenty of

water and seek medical attention. In case of contact with skin, wash immediately with plenty of soap and water.

SUPPLIED AS

Dycal in CFC-free Pressure Packs

611.07.002 Standard Package, ivory
1 x 29 g base
1 x 26 g catalyst
mixing pads

Dycal in tubes

611.05.501 Standard Package, ivory
13 g base + 11 g catalyst
mixing pads

- 611.06.501 Standard Package, dentine
13 g base + 11 g catalyst
mixing pads
- 611.05.606 6-pack, ivory
6 x 13 g base + 6 x 11 g catalyst
mixing pads
- 611.06.606 6-pack, dentine
6 x 13 g base + 6 x 11 g catalyst
mixing pads

If you have any questions, please contact:

Dentsply DeTrey GmbH
De-Trey-Strasse 1
78467 Konstanz
Germany
(0 75 31) 5 83-0
www.dentsplysirona.com

Distributor:
DENTSPLY Limited
Hamm Moor Lane
Addlestone, Weybridge
Surrey KT15 2SE, UK
Phone (0 19 32) 85 34 22

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