FluoroPost™

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FluoroPost™

Endodontic Fiber Posts, Precision Drills, Peeso Reamers

DIRECTIONS FOR USE – ENGLISH

Caution: For dental use only.

In the USA: Rx only.

1. PRODUCT DESCRIPTION

FluoroPost™ Endodontic Fiber Post is a Zirconium enriched Glass fiber/epoxy resin post designed for use in conjunction with a direct, adhesive core build-up procedure. This translucent post combines strength and elasticity and is ideal in situations where esthetics and show-through of the post may be of primary concern.

1.1 Delivery forms

FluoroPost™ Fiber Post is available in four sizes, with corresponding sized preparation instruments, color coded for easy identification.

Color Code	No. 1 Yellow	No. 2 Red	No. 3 Blue	No. 4 Green
Peeso Reamer	No. 1	No. 1	No. 2	No. 2
Precision Drill	No. 1	No. 2	No. 3	No. 4
Endo post Tip Diameter (mm)	0.80	0.80	1.00	1.00
Endo post Head diameter (mm)	1.35	1.47	1.67	1.83

1.2 Composition

FluoroPost™ Fiber Post

Zirconium enriched Glass fiber 64% volume; Epoxy Resin 36% volume

FluoroPost™ Peeso Reamers and Precision Drills
Stainless Steel

1.3 Indications

FluoroPost™ Fiber Post is indicated to augment retention of the coronal restoration in case of insufficient residual tooth structure (<4mm).

1.4 Contraindication

FluoroPost™ Fiber Post is contraindicated for use when insufficient coronal dentin remains. At least 2mm of tooth structure is required around the preparation.

1.5 Compatible cements and restoratives

• FluoroPost™ Fiber Post is to be used in conjunction with suitable adhesive cements and core build-up materials.
• FluoroPost™ Fiber Post is compatible with numerous dentin/enamel adhesive systems, self/dual cured cements and composite core build-up materials including all DENTSPLY adhesives designed for use with dual/self cured composite restoratives (see complete Directions for Use of selected adhesive). Other dentin/enamel adhesive systems, cements and core build-up materials may be used at the discretion and sole responsibility of the dental practitioner.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

	Safety alert symbol This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.
---	---

2.1 Warnings

- FluoroPost™ Peeso Reamers and Precision Drills are designed for use without water. As with any rotary instrument, heat will build up with prolonged contact. Use at slow speed (800-1200 rpm) with light intermittent pressure. Excessive speed may result in heat buildup and may endanger the stability of the instruments.
- Care should be taken to have good control on any rotary instrument to protect patients from injury. Contact with soft tissues (skin, gingiva, mucosa) may result in injuries. If injury occurs, ensure appropriate medical care.
- Placement of rubber dam or other suitable isolation is advised prior to operative procedures. Dental professionals and patients should wear protective eyewear while using the FluoroPost™ Peeso Reamers and Drills.

2.2 Precautions

- This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
- Reusable FluoroPost™ Peeso Reamers and Precision Drills must be disinfected, cleaned, dried and sterilized prior first use and after every use. The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Use of rubber dam is recommended.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- The FluoroPost™ Fiber Post must be cleaned with alcohol before insertion in the canal. Avoid touching the post with fingers after cleaning.
- Shortening of the post should be done outside the mouth.
- When inserting the FluoroPost™ Peeso Reamer and Precision Drill into the handpiece, insert the instruments as deeply as possible without excessive use of force into the chuck and check to ensure that they are secure.
- The handpiece must be in proper working condition.
- Excessive speeds while using instruments endanger the stability of the instruments.
- Avoid excessive lateral pressure, or overloading the instruments by using too much contact pressure or angle (increased danger of breakage or perforation).
- Immediately remove any instruments that are damaged, bent or no longer run concentrically.

Interactions

None Known.

2.4 Adverse reactions

Product misuse may cause failure resulting in respiratory, digestive, eye, skin, or mucosa injury. (See Contraindications, Warnings and Precautions)

2.3 Storage conditions

Keep out of direct sunlight and protect from moisture. Store at temperatures between 10°C/50°F and 24°C/75°F. Allow material to reach room temperature prior to use.

3. STEP-BY-STEP INSTRUCTIONS

- Pre-operative sterilization
FluoroPost™ Fiber Post should be steam autoclaved one time (134°C, 18 minutes minimum holding time) prior to clinical use. FluoroPost™ Peeso Reamers and Precision Drills are not supplied sterile, and should be steam autoclaved prior to first and each subsequent use (see Hygiene section).
- Place rubber dam or other suitable isolation technique.
- Root canal preparation.
Ensure that sufficient space exists to place the post. Remove residual gutta percha out of the pulp chamber and canal orifice.
- Select the correct FluoroPost™ Fiber Post size, according to the anatomical situation of tooth using the radiograph and the size table above. Select the FluoroPost™ Peeso Reamer and Precision Drill corresponding to the selected FluoroPost™ Fiber Post size.

	Rotary Instruments <ul style="list-style-type: none">Be sure instrument is securely seated in the handpiece before introducing into the patient field. Do not exert excessive pressure to prevent instrument damage or unintended removal of tooth structure/root perforation.
---	--

- Determine the length of the coronal build-up which will account for 1/3 of the final post length. Remove the root canal obturation with the selected FluoroPost™ Peeso Reamer (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted at 2/3 of its length. At least 4mm of root canal filling should remain in the apical region. In the case of curved canals, the length has to be reduced.
- Drill the canal with the selected FluoroPost™ Precision Drill (rotation speed 1000-1200 rpm).
- Trial fit the fiber post and check that it is properly seated in the canal.
- Shorten the fiber post (at the coronal end) to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- Clean the fiber post with alcohol.
- Condition and/or apply adhesive to the root canal dentin per manufacturer’s Directions for Use for the selected adhesive and/or cement.
- Apply adhesive and/or cement to the FluoroPost™ Fiber Post and post space per manufacturer’s Directions for Use for the selected adhesive and/or cement.
- Place the FluoroPost™ Fiber Post into final position and stabilize.
- Proceed with core build-up as directed by cement and/or core build-up manufacturer’s Directions for Use.

4. HIGIENE

Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

4.1 Disinfection, cleaning and sterilization

For hygiene and cross contamination caution, all instruments not marked “sterile” must be cleaned, disinfected and sterilized before each usage to prevent contamination according to ISO 17664. This concerns the first use as well as the subsequent ones.

General recommendations

- Single use marked instruments are not approved for re-use.
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments.
- For your own safety, please wear personal protective equipment (gloves, glasses).
- Use only disinfecting solution which is approved for their efficacy (DGHM-listing, CE marking, FDA approval).
- After 5 cycles, degradation of marking may occur.

Step-by-step procedure
▲ Instruments: First Usage
■ Instruments: Following Usage
✳ FluoroPost™ Fiber Post

<p>Operation</p> <p>1. Pre-Disinfection or Decontamination</p> <p>■</p>	<p>Operating mode</p> <ul style="list-style-type: none">Soak immediately just after usage all instruments in a disinfectant solution combined with proteolytic enzyme if possible	<p>Warning</p> <ul style="list-style-type: none">Follow instructions and observe concentrations and immersion time given by the manufacturer (an excessive concentration may cause corrosion or other defects on instruments). The disinfectant solution should be aldehyde free (to avoid blood impurities fixation). Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (see General Recommendation). For visible impurities are observed on instruments a pre-cleaning is recommended by brushing them manually with soft material.
<p>2a. Automated Cleaning/ Disinfection</p> <p>▲ ■ ✳</p>	<ul style="list-style-type: none">Disassemble the devices (silicone stops to be removed). Place them in a kit, support or container. Put them in the thermodisinfectior (for at least 10 min at 93°C or Ao value >3000).	<ul style="list-style-type: none">Discard any instruments with large obvious defects (broken, bent). Avoid any contact between instruments or posts when placing in the thermodisinfectior. Follow instructions and observe concentrations given by the manufacturer (see also General recommendation). The thermodisinfectior is not recommended for instruments made of aluminum, tungsten carbide or carbon steel. In case of chemical disinfection, there is a danger of remnants of the disinfectant on the instruments. Sufficient rinsing step should be available in the program. Use purified water for rinsing (max 10 germs/mL and max 0.25 endotoxin units/mL) like Aqua Purificata. Use filtered air for drying the instruments or posts. Use only approved thermodisinfectior according to EN ISO 15883, maintain and calibrate it regularly. If possible, an automated procedure if preferred.
<p>2b. Manual Cleaning/ Disinfection</p> <p>▲ ■ ✳</p>	<ul style="list-style-type: none">Disassemble the devices (silicone stops to be removed). Immerse them in the disinfectant solution assisted with ultrasonic device if suitable Rinse the devices thoroughly with clean, demineralized or distilled water and then dry them with filtered or compressed air.	<ul style="list-style-type: none">No visible impurities should be observed on the instruments. Discard any instruments with large obvious defects (broken, bent, and twisted). Avoid any contact between instruments or posts when placing in the solution, use kits, supports or container. Follow instructions and observe concentrations and time given by the manufacturer (see also General recommendation). Use purified water for rinsing (max 10 germs/mL and max 0.25 endotoxin units/mL) like Aqua Purificata. If a disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before autoclaving.
<p>3. Inspection</p> <p>■ ✳</p>	<ul style="list-style-type: none">Inspect devices and sort out those with defects. Assemble the devices (stops).	<ul style="list-style-type: none">Dirty instruments must be cleaned and disinfected again. Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of color coding or marking) affecting the resistance, safety or the performance of the instrument or posts. Protect carbon steel bur with corrosion inhibitor before packaging.
<p>4. Packaging</p> <p>▲ ■ ✳</p>	<ul style="list-style-type: none">Pack the devices in “Sterilization pouches”	<ul style="list-style-type: none">Check the validity period of the pouch given by the manufacturer to determine the shelf life. Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607.
<p>5. Sterilization</p> <p>▲ ■ ✳</p>	<ul style="list-style-type: none">Steam sterilize at: 134°C (237°F) during 18 min. for endodontic instruments and root canal filling instruments. 134°C (237°F) during at least 5 min. for all other devices	<ul style="list-style-type: none">The instruments, posts and the plastic supports must be sterilized according to the packaging labeling. Use fractionated vacuum or gravity (less preferred) autoclaves (according to EN 13060, EN 285). Use validated sterilization procedure according to ISO 17665-1. Respect maintenance procedure of the autoclave device given by the manufacturer. Use only the listed sterilization procedures.
<p>6. Storage</p> <p>▲ ■ ✳</p>	<ul style="list-style-type: none">Keep devices in sterilization packaging in a dry and clean environment.	<ul style="list-style-type: none">Sterility cannot be guaranteed if packaging is open, damaged or wet (check the packaging before using the instruments).

5. LOT NUMBER

The following numbers should be quoted in all correspondences:

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All rights reserved.