Peripac[®]

Periodontal Dressing

Parodontalverband

Pansement parodontal

Impacco parodontale

Apósito periodontal

Instructions for Use_____	English	2
Gebrauchsanweisung _____	Deutsch	6
Mode d'emploi _____	Français	10
Istruzioni per l'uso_____	Italiano	14
Instrucciones de uso _____	Español	18

Peripac®

Periodontal Dressing

CAUTION: For dental use only.

USA: Rx only

Content	Page
1 Product description	2
2 Safety notes	3
3 Step-by-step instructions	4
4 Hygiene	5
5 Lot number and expiration date	5

1 Product description

Peripac® periodontal dressing is a well proven concept in the protection of periodontal surgical sites. This single component, moisture activated periodontal dressing eliminates time-consuming mixing of paste. The dressing is peach colored and has a banana aroma.

Peripac® periodontal dressing is based on calciumsulfate-polymethylmethacrylate.

Peripac® periodontal dressing is a putty like paste which sets after water uptake from saliva.

1.1 Indications

Periodontal dressing after:

- Re-attachment (flap) surgery.
- Gingival reposition.
- Deep curettage.
- Gingivectomies and papillectomies.

Periodontal dressing is also a valuable help for:

- Protection of non-specific lesions or sutured margins.
- Temporary re-basing of immediate dentures in periodontal surgery.

1.2 Contraindications

- Use with patients who have a history of severe allergic reaction to any of the components.

1.3 Delivery forms

- Peripac® is available in 40 g jar.

1.4 Composition

- Calcium sulfate hemihydrate
- Zinc sulfate
- Zinc oxide
- Polymethylmethacrylate
- Dimethoxytetraethylenglycol
- Ascorbic acid
- Aroma
- Iron oxide pigments

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

- None known.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective measures for the dental team and patients such as glasses in accordance with local best practice.

2.3 Adverse reactions

- None known.

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store in a tightly closed container at temperatures between 10 °C and 24 °C (50 °F and 75 °F). Use product at room temperature.
- Store in a dry environment (< 80% relative humidity).
- Storage of the material under refrigerated conditions may lead to water uptake thus causing premature partial setting of the paste.
- Storage at higher temperatures will increase viscosity and thus stickiness of the material.
- Do not use after expiration date.

3. Step-by-step instructions

3.1 Preparation

1. Protect acrylic crowns with a film of Vaseline®¹ to prevent sticking.
2. Take a small quantity of the material from the jar with a dry, sterile spatula. Deposit the material on a paper napkin.
3. Tightly close jar immediately after use.

3.2 Forming the dressing

Premature hardening of the material.

The material starts to harden as soon as it comes into contact with water.

1. Avoid humidity.
2. Only form the dressing with dry fingers.

1. Form the material into a small roll of the length and thickness of the gum section to be covered.

3.3 Application

Material extending the wound or/and the periodontal margin area.

Mechanical irritation.

1. Only cover the wound area or/and the periodontal margin area with the dressing.
2. Do not cover the mobile alveolar mucosa.
3. In areas exposed to heavy masticatory function, secure the edge of the dressing lying on the crown of the tooth by spot application of a tissue adhesive.

NOTICE

Stickiness of dressing.

Dressing sticks to hand instruments.

1. Do not pack the dressing directly with an instrument.

When covering long and wide wounds or sutured margins, the roll of material is best applied with fingers.

1. Press the material well into the interdental spaces, working from the apical region towards the occlusal or incisal edges of the teeth. Technique Tip: Moisten the glove well with water or saliva will prevent the material sticking to them. Use a small, slightly moistened cotton-wool pellet held with tweezers for final adaptation into the interdental spaces.
2. Application of the dressing should not take more than 2 to 3 minutes, as further molding or addition of material is impossible once hardening has started.
3. Remove surplus material from the occlusal surface.
4. Before final hardening of the material, check the occlusion.

¹ Not a registered trademark of Dentsply Sirona, Inc.

Hardening

Premature hardening of the material.

The material starts to harden as soon as it comes into contact with water.

1. Avoid humidity.
2. Only form the dressing with dry fingers.

The hardening of the material is completed after 20 minutes.

After-treatment

A correctly applied dressing can be left without changing for 8 to 10 days if the mouth is rinsed daily with a chlorhexidine preparation. It is recommended to renew large dressing after 3 to 5 days.

1. Inform the patient about the chlorhexidine preparation.

4 Hygiene

Cross-contamination.

Infection.

1. In order to prevent cross-contamination do not touch storage container, and Instructions for Use with contaminated hands. Protect from aerosols.

5 Lot number () and expiration date ()

1. Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date

© Dentsply Sirona 2018-04-13

[These Instructions for Use are based on Master Version 03]

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¹ No es una marca comercial registrada de Dentsply Sirona, Inc.
Dentsply DeTrey GmbH