Surefil one™

Self-Adhesive Composite Hybrid

Selbstadhäsives Komposithybrid

Composite hybride auto-adhésif

Composito ibrido autoadesivo

Composite Híbrido Autoadhesivo

Capsule Extruder 2

Instructions for Use	English	2
Gebrauchsanweisung	Deutsch	10
Mode d'emploi	Français	18
Istruzioni per l'uso	Italiano	26
Instrucciones de uso	Español	34

Surefil one™

Self-Adhesive Composite Hybrid

Capsule Extruder 2

CAUTION: For dental use only.

USA: Rx only

Content	Page
1 Product description	2
2 Safety notes	3
3 Step-by-step instructions	4
4 Hygiene	7
5 Lot number, expiration date and correspondence	9

1 Product description

Surefil one™ is a self-adhesive composite hybrid offering the following features:

- Self-adhesiveness to enamel and dentin
- Fluoride release
- Bulk fill and dual cure

It's radiopacity is equivalent to 2 mm aluminum. (The radiopacity of aluminum is comparable to dentin.)

Surefil one™ meets the requirements regarding flexural strength for direct composite restorations¹.

Capsule Extruder 2 is the device designed for extrusion and direct application of Surefil one™ restorative capsules.

1.1 Indications

- Class I to Class V restorations.

1.2 Contraindications

- Direct or indirect pulp-capping.
- Use in patients with known allergies to any of the components.

1.3 Delivery forms

- Pre-dosed mixing capsule for direct intraoral application with a minimal dispensable amount of 0.3 g.

¹ Based on ISO standard 4049:2009 (Type 1, class III).

1.4 Composition

- Aluminum-phosphor-strontium-sodium-fluoro-silicate glass
- Water
- Highly dispersed silicon dioxide
- Acrylic acid
- Polycarboxylic acid
- Ytterbium fluoride
- Bifunctional acrylate
- Self cure Initiator
- Iron oxide pigments
- Barium sulfate pigment
- Manganese pigment
- Camphorquinone (photoinitiator)
- Stabilizer

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

The mixed material contains polycarboxylic acid and may be irritating to skin, eyes and oral mucosa and may cause sensitization by skin contact and allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, remove material with cotton and wash thoroughly with soap and water.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use. Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective eyewear, mask, clothing and gloves for the dental team and protective eyewear for patients.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Capsules are intended for single use only. Discard after use. Do not reuse.
- If material has started to set excessive pressure may result in rupture of the capsule.
- Increased temperature and/or humidity reduce available working time and accelerate setting time. Material will set faster in the mouth than at the bench. Do not use bench cure rates to judge mouth cure rates.

- Interactions:
 - Do not use eugenol- or hydrogen peroxide-containing materials in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
 - The restorative is a light cured material. Therefore, it should be protected from ambient light. Proceed immediately once material has been placed.
 - Using coated matrix bands may lead to a whitish appearance of the material surface that should be finished away.

2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Sealed packs can be stored at temperatures between 2 °C and 24 °C (36 °F and 75 °F). Refrigerated storage is not recommended. Do not freeze!
- Use capsules at room temperature.
- Protect from moisture.
- Do not use after expiration date.

3 Step-by-step instructions

3.1 Shade selection

A1, A2, A3, A3.5 and Bleach White.

3.2 Cavity preparation

1. In all cavity classes cavity preparation may be kept to the minimum required for caries removal.
2. Use fine diamonds or carbides to finish cavity preparation.
3. Prepare butt joint – avoid beveled margins.
4. Rinse surface with water spray and carefully remove excess rinsing water with air spray. Do not desiccate the tooth substance – cavity surface should be moist.
5. Use cotton rolls or rubber dam to isolate the cavity from contamination.

3.3 Placement of matrix

1. In Class II cavities place an uncoated matrix (e.g., AutoMatrix® matrix system or Palodent® matrix system; see complete Instructions for Use) and wedge.
2. Use a dead soft, thin matrix band.
3. In Class V cavities cervical matrices may be used.

3.4 Pulp protection

1. In deep cavities cover the dentin close to the pulp (less than 1 mm) with a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see complete Instructions for Use) leaving the rest of the cavity surface free for bonding.

3.5 Application

Time after activation of capsule

According to ISO 9917-2²

Working time (at room temperature)	2 minutes
Setting time (at 37 °C/99 °F)	3 minutes 30 seconds

Typical clinical protocol

Mixing (capsule)	10 seconds
Working time (at 30 °C/86 °F)	1 minute 30 seconds
Light curing ³ of surface layer	20 sec ≥ 800 mW/cm ²
Setting time at 30 °C/86 °F	6 minutes (after activation)

Capsule

3.5.1 Activation of the capsule

1. Activate the capsule by pressing the capsule onto a stable surface and depressing the plunger to its limit (at this position the plunger will overlap by about 2 mm).

3.5.2 Mixing

1. **Immediately** place the activated capsule in a capsule mixer (4,200-5,000 oscillations/minute) and mix for **10 seconds unless stated otherwise in Table 1:**

Table 1: Recommended mixing times with other than standard capsule mixers (4,200-5,000 oscillations/minute)

Trademark/Manufacturer	Min (sec)	Max (sec)
Rotomix⁴ – with centrifugation 3M Oral Care	8	8

3.5.3 Application of restorative

Surefil one™ was designed to provide tooth colored restorations in posterior teeth. However, anterior restorations (class III-V) are indicated as well but shade range, opacity, and the need of butt joint margins may limit the esthetic outcome of restorations in the smile line.

1. Immediately remove the capsule from the capsule mixer and place into the Capsule Extruder 2. The capsule may be rotated 360° to gain the proper angle of entrance into the cavity. Do not apply excessive force.

² Surefil one™ is a Restorative/Cement class III, acc. to ISO 9917-2.

³ Curing light designed to cure materials containing camphorquinone (CQ) initiator. Peak of spectrum in the range of 440-480 nm.

⁴ No product of Dentsply Sirona, Inc.

2. Click the trigger of the Capsule Extruder 2 until paste is seen through the nozzle. **Immediately start dispensing.** Begin dispensing at the deepest part of the cavity, keeping tip close to cavity floor. While moving the capsule tip in all areas of the cavity, gradually withdraw tip as cavity is filled. Avoid lifting the tip out of dispensed material while dispensing to minimize air entrapment. At the completion of dispensing, wipe tip against cavity wall while withdrawing from the operative field. It is recommended to slightly overfill the cavity.
3. Start excess removal and contouring **immediately after placement.**
4. Do not manipulate the material after start of set (see Working time 3.5.4 below).
5. After application, press the release button of the Capsule Extruder 2 until the push rod releases. Slightly rotating the capsule, the empty capsule can easily be removed.

Using two capsules for large cavities

- Either apply second capsule within working time (see 3.5.4) of first capsule by placing horizontal increments so that always the freshly applied material from the second capsule is the upper most layer to allow manipulation before end of respective working time.
- Alternatively, light cure surface of material applied from first capsule before applying second capsule. This technique allows vertical layering (e.g., fill approximal box to marginal ridge first, light cure and then fill remaining class I cavity with second capsule).

3.5.4 Working and setting time

- Minimum working time in the oral cavity is **1 minute 30 seconds** from start of mixing. Operatory light or LED light from loupes will shorten working time.
- Light cure the surface layer for 20 seconds⁵ for immediate start of contouring and finishing. Otherwise, wait at least **6 minutes** (after activation) until further manipulation.

Gross contamination during placement (cavity full of rinsing water) and setting (restoration flooded by saliva).

Dehydrated or diluted self-adhesive restorative.

1. Do not use high-volume suction close to restoration to prevent desiccation.
2. Prevent saliva contamination during application and early setting.

3.5.5 Finishing

1. Keep the matrix in place for at least **6 minutes after activation.**
2. Start finishing immediately after surface light curing or no earlier than 6 minutes after activation with cups, discs, and points (Enhance[®] Finishing System).
3. Use water spray during contouring, finishing, and polishing to avoid desiccation of the surface layer.

⁵ Curing light designed to cure materials containing camphorquinone (CQ) initiator. Peak of spectrum in the range of 440-480 nm.

4 Hygiene

Cross-contamination.

Infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Reprocess reusable products according as described below.

4.1 Capsule Extruder 2

NOTICE

Wrong cleaning or disinfection method.

Damage to the Capsule Extruder 2.

1. Do not submerge into disinfection solutions.
2. Do not ultrasonic clean.
3. Start reprocessing within 1 hour after use.

Instructions for cleaning and disinfecting Capsule Extruder 2	
Warnings	<ul style="list-style-type: none">• Intermediate-level disinfection is appropriate for the Capsule Extruder 2.• Do not submerge into disinfection solutions.• Do not ultrasonic clean.• Do not allow the Capsule Extruder 2 to exceed 137 °C/279 °F.
Limitations on Reprocessing	<ul style="list-style-type: none">• Repeated reprocessing has minimum effect on these instruments. End of life is normally determined by wear and damage due to use.• The Capsule Extruder 2 can be subjected to the manual cleaning process up to 1,000 times.• The Capsule Extruder 2 can be subjected to the automated washer/disinfector process up to 200 times.• The Capsule Extruder 2 can be subjected to the steam autoclave sterilization process up to 1,000 times.• Cold liquid immersion disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.
Point of Use	<ol style="list-style-type: none">1. Right after use, remove capsule from the Capsule Extruder 2 and discard capsule according to local regulations.2. Clean Capsule Extruder 2 using a soft brush and running cold tap water.3. Discard gloves according to local regulations and disinfect hands with an appropriate bactericidal, virucidal, and fungicidal hand disinfectant solution according to local regulations. Use according to disinfectant solution manufacturer's Instructions for Use.4. Start reprocessing within 1 hour after use.
Containment and Transportation	No particular requirements. Handle with care.

Cleaning and disinfection: Automated	<p>As an alternative to manual cleaning and disinfection, the device can also be processed in a washer-disinfector. Use only properly maintained, calibrated, and approved washer-disinfector according to ISO 15883-1.</p> <ol style="list-style-type: none"> 1. Right after use, remove capsule from the device and discard capsule according to local regulations. 2. Clean Capsule Extruder 2 with a soft brush and running cold tap water. 3. Place device in the washer-disinfector allowing water and detergent to enter into and drain off from device openings during washer-disinfection program. Follow washer-disinfector's Instructions for Use. 4. Run washer-disinfection program with A0 value $\geq 3,000$ (e.g., 5 min at ≥ 90 °C/194 °F), using appropriate detergents.
Cleaning and Disinfection: Manual	<ol style="list-style-type: none"> 1. Use a new clean pair of examination gloves. 2. Clean Capsule Extruder 2 with a soft brush, tap water, and liquid soap until free of visible contamination. Pay special attention to device seams and insertions. 3. Remove cleaning solution residue with a cloth soaked with tap water. 4. Dry Capsule Extruder 2 with a lint-free, single-use cloth. 5. Manual disinfection is not recommended.
Packaging	<p>No particular requirements.</p>
Sterilization/ autoclaving	<p>After cleaning, the Capsule Extruder 2 must also be autoclaved.</p> <ol style="list-style-type: none"> 1. Steam autoclave Capsule Extruder 2 unwrapped at 134 °C/2 bar = 273 °F/29 psi for 3 minutes 30 seconds.
Drying	<ul style="list-style-type: none"> • Dry Capsule Extruder 2 with a lint-free single-use cloth. • Allow the Capsule Extruder 2 to fully air dry before storage.
Maintenance	<ul style="list-style-type: none"> • Visually inspect to ensure that all contamination has been removed. • Regularly apply dental instrument lubricant. • Do not continue to use damaged or soiled Capsule Extruder 2.
Storage	<ul style="list-style-type: none"> • Store the Capsule Extruder 2 at room temperature, away from moisture or excessive humidity.
Manufacturer Contact	<p>Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.</p>

The instructions provided above have been validated by Dentsply Sirona as being capable of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing is actually performed using equipment, materials and trained personnel in the processing facility to achieve the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Use of other reprocessing procedures/methods is at the discretion and sole responsibility of the practitioner.

Reprocessing Validation Conditions		
Component	Process	Method
Capsule Extruder 2	Automated washer/ disinfecter	Miele Washer-Disinfecter G7835 CD using program vario TD. Washer-disinfection program with A0 value $\geq 3,000$ (e.g., 5 min at ≥ 90 °C/194 °F). Detergents used for validation: neodisher® MediClean [0.5%] (alkaline detergent) and neodisher® Z [0.1%] (acid neutralization and cleaning detergent); both detergents from Dr. Weigert, Hamburg, Germany. ⁶
Capsule Extruder 2	Manual Cleaning	Tap water, liquid soap.
Capsule Extruder 2	Autoclave Sterilization	Unwrapped at 134 °C/2 bar = 273 °F/29 psi for 3 minutes 30 seconds.

4.2 Disposal

Dispose of in accordance with local regulations.

5 Lot number (), expiration date () and correspondence

- Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
- The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date
- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

© Dentsply Sirona 2019-06-07

[These Instructions for Use are based on Master Version 02]

⁶ No products or registered trademarks of Dentsply Sirona, Inc.

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⁶ No son productos ni marcas comerciales registradas de Dentsply Sirona, Inc.
Dentsply DeTrey GmbH
Distributed in U.S.A. by
Dentsply LLC