# X-Post<sup>®</sup>

## Endodontic Fiber Post

Endodontischer Glasfaserstift  
Tenon endodontique fibré  
Perno endodontico in fibra  
Postes de fibra de endodoncia

Instructions for Use \_\_\_\_\_ English 2

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# X-Post®

## Endodontic Fiber Post

**Caution:** For dental use only. **In the USA:** RX only.

### Composition

- Quartz fibers 60% volume
- Epoxy resin 40% volume

### Indications for use

These instruments are to be used only in a clinical or hospital environment, by qualified users.

### Application field

Endodontic fiber post for direct, adhesive core-build up. In case of insufficient residual tooth substance (< 4 mm), the **X-Post®** may be used to support the coronal restoration.

### Contraindications

Insufficient residual dentin: At least 2 mm of tooth structure is required around the preparation.

### Warnings

None known.

### Precautions

The post must be cleaned with alcohol before insertion in the canal.

Avoid touching the posts with your fingers after cleaning.

Shortening the post should be done outside of the mouth.

Dentsply Sirona recommends the use of a rubber dam.

### Adverse Reactions

In the present technical state, no adverse reaction has been reported so far.

### Step-by-step Instructions

1. Root canal preparation. Ensure that sufficient space exists to place the post. Remove residual gutta percha out of the root canal.
2. Select correct **X-Post®** size, according to anatomical situation of tooth using the radiograph and the information provided below.  
Select the Largo® Peeso Reamer and Easypost™ Precision Drill corresponding to the selected **X-Post®** size.

X-Post®	No. 1	No. 2	No. 3	No. 4
Color Code	Yellow	Red	Blue	Green
Largo® Peeso Reamer	No. 1	No. 1	No. 2	No. 2
Easypost™ Precision Drill	No. 1	No. 2	No. 3	No. 4
Tip diameter (mm)	0.80	0.80	0.80	0.80
Head diameter (mm)	1.35	1.47	1.67	1.83

3. Determine the length of the coronal build-up which will account for 1/3 of the final post length. Remove the root canal obturation with the selected Largo® Peeso Reamer (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted at 2/3 of its length. At least 4 mm of root canal filling should remain in the apical region.  
In case of curved canals, the length has to be reduced.
4. Drill the canal with the selected Easypost™ Precision Drill (rotation speed 1000-1200 rpm).
5. Check that the post is properly seated in the canal.
6. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
7. Clean the post with alcohol.
8. Apply etchant in the canal and to the exposed dentin for 15 seconds (e.g. DeTrey® Conditioner 36). Rinse for 10 seconds. Dry gently with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.  
**Note:** phosphoric etching is optional for Prime&Bond active™, Prime&Bond universal™ and Prime&Bond elect®.<sup>1</sup>
9. Mix Prime&Bond® adhesives<sup>2</sup> and Self Cure Activator. Apply 1 coat of the mixture with a brush (e.g. Applicator Tips) in the root canal and leave for 20 seconds. Remove excess using paper points.
10. Gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat the application.
11. Light cure the adhesive mixture in the root canal for 10 seconds<sup>3</sup>.  
**Note:** Light cure step is not required when using core-X® flow + Self Cure Activator with either Prime&Bond active™, Prime&Bond universal™, Prime&Bond elect® or Prime&Bond® XP.

12. Apply a single coat of Prime&Bond® adhesive<sup>2</sup> and Self Cure Activator mixture to the post. Gently air dry for 5 seconds.

For further treatment procedure, please consult Illustrated Technique Guide and/or the Instructions for Use of core-X® flow or Calibra®.

### Disinfection, Cleaning and Sterilization

For hygiene and cross contamination caution, all instruments not marked "sterile" must be cleaned, disinfected and sterilized before each usage to prevent any contamination according to ISO 17664. This concerns the first use as well as the subsequent ones.

### General recommendation

- Single use marked instruments are not approved for re-use.
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments.
- For your own safety, please wear personal protective equipments (gloves, glasses).
- Use only disinfecting solution which is approved for their efficacy (DGHM-listing, CE marking, FDA approval).
- Tungsten Carbide burs, Plastic Support, Hand instruments and NiTi instruments are degraded by Peroxide hydrogen (H<sub>2</sub>O<sub>2</sub>) solution.
- NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5%.
- Instruments made of aluminium are degraded in presence of caustic soda solutions with mercury salt. Do not use acid (pH < 6) or alkaline (pH > 8) solutions.
- After 5 cycles, degradation of marking may appear.

<sup>1</sup> Some adhesives are not available in all countries.

<sup>2</sup> Prime&Bond active™, Prime&Bond universal™, Prime&Bond elect®, Prime&Bond® XP (in some countries sold as XP Bond®), prime&bond® one Etch&Rinse, Prime&Bond® NT.

<sup>3</sup> Check curing light for minimum light output of at least 800 mW/cm<sup>2</sup>. Cure for at least 20 seconds if light output is between 500 and 800 mW/cm<sup>2</sup>.

Step-by-step procedure

Implantable Devices

Instruments: Following usage

Instruments: First usage

A. Devices

	Operation	Operating mode	Warning				
1.	<b>Pre-Disinfection or Decontamination</b>	- Soak immediately just after usage all instruments in a disinfectant solution combined with proteolytic enzyme if possible.	- Follow instructions and observe concentrations and immersion time given by the manufacturer (an excessive concentration may cause corrosion or other defects on instruments). - The disinfectant solution should be aldehyde free (to avoid blood impurities fixation). - Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (see <b>General recommendation</b> ). - For visible impurities are observed on instruments a pre-cleaning is recommended by brushing them manually with soft material.			X	
2a.	<b>Automated Cleaning / Disinfection</b>	- Disassemble the devices (silicone stops to be removed). - Place them in a kit, support or container. - Put them in the thermodisinfectant (for at least 10 min at 93 °C or Ao value > 3000).	- Discard any instruments with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the thermodisinfectant. - Follow instructions and observe concentrations given by the manufacturer (see also <b>General recommendation</b> ). - The thermodisinfectant is not recommended for instruments made of aluminium, tungsten carbide or carbon steel. In case of chemical disinfection, there is a danger of remnants of the disinfectant on the instruments. - Sufficient rinsing step should be available in the program. Use purified water for rinsing (max 10 germs/ml and max 0.25 endotoxin units/ml) like Aqua Purificata. - Use filtered air for drying the instruments or posts. - Use only approved thermodisinfectant according to EN ISO 15883, maintain and calibrate it regularly. - If possible, an automated procedure is preferred.		X	X	X
2b.	<b>Manual Cleaning / Disinfection</b>	- Disassemble the devices (silicone stops to be removed). - Immerse them in the disinfectant solution assisted with ultrasonic device if suitable. - Rinse the devices thoroughly with clean, demineralised or distilled water and then dry them with filtered compressed air.	- No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Avoid any contact between instruments or posts when placing in the solution, use kits, supports or container. - Follow instructions and observe concentrations and time given by the manufacturer (see also <b>General recommendation</b> ). - Use purified water for rinsing (max 10 germs/ml and max 0.25 endotoxin units/ml) like Aqua Purificata. - If a disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before autoclaving.		X	X	X
3.	<b>Inspection</b>	- Inspect devices and sort out those with defects. - Assemble the devices (stops).	- Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the instrument or posts. - Protect carbon steel bur with corrosion inhibitor before packaging.			X	X
4.	<b>Packaging</b>	- Pack the devices in "Sterilisation pouches".	- Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141 °C (286 °F) and in accordance with EN ISO 11607.		X	X	X
5.	<b>Sterilization</b>	- Steam sterilisation at: 134 °C (237 °F) during 18 min. for endodontic instruments and root canal filling instruments. 134 °C (237 °F) during at least 5 min. for all other devices.	- The instruments, posts and the plastic supports must be sterilized according to the packaging labelling. - Use fractionated vacuum or gravity (less preferred) autoclaves (according to EN 13060, EN 285). - Use validated sterilisation procedure according to ISO 17665-1. - Respect maintenance procedure of the autoclave device given by the manufacturer. - Use only the listed sterilization procedures.		X	X	X
6.	<b>Storage</b>	- Keep devices in sterilization packaging in a dry and clean environment.	- Sterility cannot be guaranteed if packaging is open, damaged or wet (check the packaging before using the instruments).		X	X	X

B. Filling material

	Operation	Operating mode	Warning
1.	<b>Disinfection</b>	- Immerse the obturation devices in alcohol during 30 sec at ambient temperature.	- Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (see <b>General recommendation</b> ).

© Dentsply Sirona 2017-11

[These Instructions for Use are based on Master Version "X POST RADIX FIBER POST\_DFU\_1117\_MASTER\_DSE\_EN Rev 00"]

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Dentsply DeTrey GmbH